

Single and double head injectors for computed tomography (21 bar/305 psi)



Article number

CT2/MRT Set ELS

31 7625

2 pieces	ELS 200 ml syringe	31 6026
2 pieces	QFT	31 6040
1 piece	MRS 222 MR-tube system	31 7100
1 piece	ES 224/150 patient line	31 8151

Packing unit: 20 pieces/carton – 600 pieces/pallet



CT2/MRT Set ELS (S)

31 7627

2 pieces	ELS 200 ml syringe	31 6026
1 piece	MRS 222 MR-tube system	31 7100
1 piece	ES 224/150 patient line	31 8151

Packing unit: 20 pieces/carton – 600 pieces/pallet



ELS 200 ml syringe, QFT

31 6025

Residual volume: 1.5 ml

ELS 200 ml (S) syringe

31 6026

Residual volume: 1.5 ml

Packing unit: 50 pieces/carton – 1,500 pieces/pallet



MRS 222 MR-tube system 31 7099 with inline drip chamber (for Scanbag CM-supply bag)

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.6 ml
NaCl-Side	36 cm	1.3 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet

CONSUMABLES : CT-D

Single and double head injectors for computed tomography (21 bar/305 psi)



MRS 222 MR-tube system

Article number

31 7100

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.6 ml
NaCl-Side	36 cm	1.3 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet



Double head tube system with inline drip chamber (for Scanbag CM-supply bag)

31 7109

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.4 ml
NaCl-Side	36 cm	1.1 ml
Patient line (ID: 1.5 mm)	150 cm	2.6 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet



MRS 225 MR-tube system

31 7125

Inner diameter: 2.0 mm

	Length	Filling volume
CM-Side	6 cm	0.2 ml
NaCl-Side	32 cm	1.0 ml

Packing unit: 50 pieces/carton – 4,000 pieces/pallet



ES 223 tube system with spike

31 7080

	ID	Length	Filling volume
Pressure side	2.0 mm	13 cm	0.4 ml
Suction side	2.7 mm	100 cm	6,3 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet

CONSUMMABLES : CT-D

Single and double head injectors for computed tomography (21 bar/305 psi)



Article number

EK 223 tube system with drip chamber **31 7081**

	ID	Length	Filling volume
Pressure side	2.0 mm	13 cm	0.4 ml
Suction side	2.7 mm	100 cm	5.8 ml
drip chamber			10 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet

just for
Accutron CT



ES 224/25 patient line with valve **31 8020**

Length: 25 cm
Innen diameter: 1.5 mm
Filling volume: 0.4 ml

Packing unit: 100 pieces/carton – 8,000 pieces/pallet



1-way stop cock with 30 cm extension line, luer lock **31 8025**

Length: 32 cm
Inner diameter: 2.0 mm
Filling volume: 1.2 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Patient line 150 cm with terminal check valves **31 8026**

Length: 150 cm
Inner diameter: 2.0 mm
Filling volume: 4.8 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMMABLES : CT-D

Single and double head injectors for computed tomography (21 bar/305 psi)




Patient line 150 cm with two valves

Article number

31 8036

Length: 150 + 10 cm
Inner diameter: 2.0 mm
Filling volume: 5.2 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



DW 229 connector with 3-way stop cock

31 8042

Filling volume: 0.5 ml

Packing unit: 100 pieces/carton



ES 224 Patient line

Inner diameter: 1.5 mm

ES224/80, length: 80 cm, filling volume: 1.4 ml

31 8081

ES224/100, length: 100 cm, filling volume: 1.8 ml

31 8101

ES224/150, length: 150 cm, filling volume: 2.7 ml

31 8151

ES224/200, length: 200 cm, filling volume: 3.6 ml

31 8199

ES224/250, length: 250 cm, filling volume: 4.4 ml

31 8250

ES224/300, length: 300 cm, filling volume: 5.3 ml

31 8300

Packing unit: 200 pieces/carton – 11,200 pieces/pallet



Patient line with valve

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.6 ml

31 8152

Length: 200 cm, filling volume: 3.5 ml

31 8198

Length: 250 cm, filling volume: 4.4 ml

31 8251

Length: 300 cm, filling volume: 5.3 ml

31 8301

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMMABLES : CT-D

Single and double head injectors for computed tomography (21 bar/305 psi)

Article number



Patient line with inline valve

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.6 ml

31 8154

Length: 200 cm, filling volume: 3.5 ml

31 8204

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line with 2 valves

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.6 ml

31 8217

Length: 200 cm, filling volume: 3.5 ml

31 8218

Length: 250 cm, filling volume: 4.4 ml

31 8219

Length: 300 cm, filling volume: 5.3 ml

31 8220

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line with valve

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

31 8253

Length: 300 cm, filling volume: 5.3 ml

31 8303

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

31 8252

Length: 300 cm, filling volume: 5.3 ml

31 8302

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMMABLES : CT-D

Single and double head injectors for computed tomography (21 bar/305 psi)

Article number



Coiled line with inline valve

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

Length: 300 cm, filling volume: 5.3 ml

31 8254

31 8304

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line 150 cm

31 8213

Length: 150 cm

Inner diameter: 1.6 mm

Filling volume: 3.1 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line with valve

Inner diameter: 1.6 mm

Coiled line 150 cm with valve

31 8214

Length: 150 cm, filling volume: 3.1 ml

Coiled line 180 cm with valve

31 8181

Length: 180 cm, filling volume: 3.7 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line 235 cm with Y-connector with check valve

31 8235

Length: 235 cm + 1 x 10 cm + 1 x 12 cm

Inner diameter: 1.6 mm

Total Filling volume: 5.6 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMMABLES : CT-D

Single and double head injectors for computed tomography (21 bar/305 psi)

Article number



SD 225 suction tube with 3-way stop cock and spike 31 9080

Filling volume: 0.3 ml (3-way stop cock) + 6,3 ml (suction tube)

Packing unit: 50 pieces/carton – 4,000 pieces/pallet



DW 226 3-way stop cock 31 9085

Filling volume: 0.3 ml

Packing unit: 300 pieces/carton



SP 227 vented spike, luer lock 31 9095

Filling volume: 0.4 ml

Packing unit: 50 pieces/carton



Minispike, green, luer lock 31 9099

Filling volume: 0.2 ml

Packing unit: 300 pieces/carton