

## CONSUMMABLES : HP

### Single and double head high-pressure injectors (83 bar/1,200 psi)



#### ELS 200 ml syringe, QFT

Residual volume: 1.5 ml  
At 83 bar only usable with pressure sleeve!

Article number  
31 6025

#### ELS 200 ml syringe

At 83 bar only usable with pressure sleeve!

31 6026

Packing unit: 50 pieces/carton – 1,500 pieces/pallet



#### HS 224 RA Patient line

Inner diameter: 1.8 mm  
Rotating Adapter  
Braided PU tube with check valve

HS 224/100 RA, length: 100 cm, filling volume: 2.5 ml

31 5100

HS 224/120 RA, length: 120 cm, filling volume: 3.1 ml

31 5126

HS 224/150 RA, length: 150 cm, filling volume: 3.9 ml

31 5136

HS 224/180 RA, length: 180 cm, filling volume: 4.6 ml

31 5135

Packing unit: 50 pieces/carton



#### HS 224 Patient line

Inner diameter: 1.8 mm  
Braided PU tube with check valve

HS 224/120, length: 120 cm, filling volume: 3.1 ml

31 5125

HS 224/180, length: 180 cm, filling volume: 4.6 ml

31 5134

Packing unit: 50 pieces/carton



#### HS 224 RA/180 Patient line with valve

31 5124

Inner diameter: 1.8 mm  
Rotating Adapter  
Braided PU tube with check valve with check valve

Length: 18 cm  
Filling volume: 4.6 ml

Packing unit: 50 pieces/carton

## CONSUMMABLES : HP

### Single and double head high-pressure injectors (83 bar/1,200 psi)



#### HP 180 Patient line PVC

Article number

31 5132

Inner diameter: 1.5 mm

Length: 180 cm,  
Filling volume: 3.2 ml

Packing unit: 50 pieces/carton



#### HP 150 RA Patient line PVC-flex

31 5138

Inner diameter: 1.5 mm  
Rotating Adapter

Length: 150 cm  
Filling volume: 2.7 ml

Packing unit: 50 pieces/carton



#### HP Patient line with valve, PVC-flex

Inner diameter: 1.5 mm  
Rotating Adapter/PVC-flex tube, with check valve

HP 120 RA, length: 120 cm, filling volume: 2.2 ml  
HP 150 RA, length: 150 cm, filling volume: 2.7 ml

31 5137

31 5121

Packing unit: 50 pieces/carton



#### HP 110 RA Patient line, PVC-flex

31 5139

Total length: 110 cm  
Inner diameter: 1.5 mm  
Filling volume: 2.0 ml

Packing unit: 25 pieces/carton



# CONSUMMABLES : HP

## Single and double head high-pressure injectors (83 bar/1,200 psi)



HP 110 RA Patient line with valve,  
PVC-flex

Total length: 110 cm  
Inner diameter: 1.5 mm  
Filling volume: 2.0 ml

Packing unit: 25 pieces/carton

Article number

31 5140



HSD 525 suction tube

Length: 100 cm  
Inner diameter: 2.7 mm  
Filling volume: 0.3 ml 3-way stop cock,  
6.0 ml tube + spike

Packing unit: 50 pieces/carton

31 5080



HSD 527 suction tube with valve

Length: 100 cm  
Inner diameter: 2.7 mm  
Filling volume: 0.2 ml valve,  
16 ml tube with drip chamber

Packing unit: 50 pieces/carton

31 5082



Accutron HP-D tube system

	ID	Length	Filling volume
Pressure side	1.8 mm	26.5 cm	1.4 ml (2 x 0.7 ml)
Suction side	2.7 mm	100 cm	16 ml

Packing unit: 25 pieces/carton

31 7183

## CONSUMMABLES : HP

### Single and double head high-pressure injectors (83 bar/1,200 psi)



**ScanBag® HP-D tube system**  
(for Scanbag CM-supply bag)

Article number

**31 7184**

	ID	Length	Filling volume
Pressure side	1.8 mm	26.5 cm	1.4 ml (2 x 0.7 ml)
Suction side	2.7 mm	100 cm	16 ml

Packing unit: 25 pieces/carton

just for  
Accutron HP-D



**SP 227 vented spike, luer lock**

**31 9095**

Filling volume: 0.4 ml

Packing unit: 50 pieces/carton



**HDW 526 high-pressure 3-way stop cock**

**31 5085**

Filling volume: 0.3 ml

Packing unit: 48 pieces/carton



**Minispike, green, luer lock**

**31 9099**

Filling volume: 0.2 ml

Packing unit: 300 pieces/carton