

Single and double head high-pressure injectors (83 bar/1,200 psi)



Article number

ELS 200 ml syringe, QFT

31 6025

Residual volume: 1.5 ml
At 83 bar only usable with pressure sleeve!

ELS 200 ml syringe

31 6026

At 83 bar only usable with pressure sleeve!

Packing unit: 50 pieces/carton – 1,500 pieces/pallet



HS 224 RA Patient line

Inner diameter: 1.8 mm
Rotating Adapter
Braided PU tube with check valve

HS 224/100 RA, length: 100 cm, filling volume: 2.5 ml

31 5100

HS 224/120 RA, length: 120 cm, filling volume: 3.1 ml

31 5126

HS 224/150 RA, length: 150 cm, filling volume: 3.9 ml

31 5136

HS 224/180 RA, length: 180 cm, filling volume: 4.6 ml

31 5135

Packing unit: 50 pieces/carton



HS 224 Patient line

Inner diameter: 1.8 mm
Braided PU tube with check valve

HS 224/120, length: 120 cm, filling volume: 3.1 ml

31 5125

HS 224/180, length: 180 cm, filling volume: 4.6 ml

31 5134

Packing unit: 50 pieces/carton



HS 224 RA/180 Patient line with valve

31 5124

Inner diameter: 1.8 mm
Rotating Adapter
Braided PU tube with check valve with check valve

Length: 18 cm
Filling volume: 4.6 ml

Packing unit: 50 pieces/carton

CONSUMMABLES : HP-D

Single and double head high-pressure injectors (83 bar/1,200 psi)



HP 180 Patient line PVC

Article number

31 5132

Inner diameter: 1.5 mm

Length: 180 cm,
Filling volume: 3.2 ml

Packing unit: 50 pieces/carton



HP 150 RA Patient line PVC-flex

31 5138

Inner diameter: 1.5 mm
Rotating Adapter

Length: 150 cm
Filling volume: 2.7 ml

Packing unit: 50 pieces/carton



HP Patient line with valve, PVC-flex

Inner diameter: 1.5 mm
Rotating Adapter/PVC-flex tube, with check valve

HP 120 RA, length: 120 cm, filling volume: 2.2 ml
HP 150 RA, length: 150 cm, filling volume: 2.7 ml

31 5137

31 5121

Packing unit: 50 pieces/carton



HP 110 RA Patient line, PVC-flex

31 5139

Total length: 110 cm
Inner diameter: 1.5 mm
Filling volume: 2.0 ml

Packing unit: 25 pieces/carton

straight
packed

CONSUMMABLES : HP-D

Single and double head high-pressure injectors (83 bar/1,200 psi)

Article number



HP 110 RA Patient line with valve,
PVC-flex

31 5140

Total length: 110 cm
Inner diameter: 1.5 mm
Filling volume: 2.0 ml

Packing unit: 25 pieces/carton



HSD 525 suction tube

31 5080

Length: 100 cm
Inner diameter: 2.7 mm
Filling volume: 0.3 ml 3-way stop cock,
6.0 ml tube + spike

Packing unit: 50 pieces/carton



HSD 527 suction tube with valve

31 5082

Length: 100 cm
Inner diameter: 2.7 mm
Filling volume: 0.2 ml valve,
16 ml tube with drip chamber

Packing unit: 50 pieces/carton



Accutron HP-D tube system

31 7183

	ID	Length	Filling volume
Pressure side	1.8 mm	26.5 cm	1.4 ml (2 x 0.7 ml)
Suction side	2.7 mm	100 cm	16 ml

Packing unit: 25 pieces/carton

just for
Accutron HP-D

CONSUMMABLES : HP-D

Single and double head high-pressure injectors (83 bar/1,200 psi)



ScanBag® HP-D tube system (for Scanbag CM-supply bag)

Article number

31 7184

	ID	Length	Filling volume
Pressure side	1.8 mm	26.5 cm	1.4 ml (2 x 0.7 ml)
Suction side	2.7 mm	100 cm	16 ml

Packing unit: 25 pieces/carton

just for
Accutron HP-D



SP 227 vented spike, luer lock

31 9095

Filling volume: 0.4 ml

Packing unit: 50 pieces/carton



HDW 526 high-pressure 3-way stop cock

31 5085

Filling volume: 0.3 ml

Packing unit: 48 pieces/carton



Minispike, green, luer lock

31 9099

Filling volume: 0.2 ml

Packing unit: 300 pieces/carton