

Double head injectors for magnetic resonance imaging (21 bar/305 psi)

Article number



MR-Set ELS 65/65 ml

31 7602

2 pieces	ELS 65 ml syringe	31 6065
1 piece	MR-tube system with valves	31 7103
1 piece	standard spike	
1 piece	minispike	

Packing unit: 50 pieces/carton

Suitable for 65/65 ml version



MR-Set ELS 65/200 ml

31 7603

1 piece	ELS 65 ml syringe	31 6065
1 piece	ELS 200 ml syringe	31 6025
1 piece	Double head tube system	31 7101
1 piece	SP 227 vented spike, luer lock	31 9095
1 piece	Minispike, green, luer lock	31 9099

Packing unit: 24 pieces/carton – 384 pieces/pallet

Suitable for 65/200 ml version



MR-Kit ELS 65 ml

31 7604

1 piece	ELS 65 ml syringe	31 6065
1 piece	MR-tube system with valves	31 7103
1 piece	standard spike	

Packing unit: 50 pieces/carton



CT2/MRT Set ELS

31 7625

2 pieces	ELS 200 ml syringe	31 6026
2 pieces	QFT	31 6040
1 piece	MRS 222 MR-tube system	31 7100
1 piece	ES 224/150 patient line	31 8151

Packing unit: 20 pieces/carton – 600 pieces/pallet

Suitable for 200/200 ml version

just for
Accutron MR

CONSUMMABLES : MR

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



CT2/MRT Set ELS (S)

Article number

31 7627

2 pieces ELS 200 ml syringe
1 piece MR 222 MR-tube system
1 piece ES 224/150 patient line

31 6026
31 7100
31 8151

Packing unit: 20 pieces/carton – 600 pieces/pallet

Suitable for 200/200 ml version

just for
Accutron MR



ELS 65 ml syringe

31 6065

Residual volume: 1.0 ml

Packing unit: 50 pieces/carton – 2,100 pieces/pallet



ELS 200 ml syringe, QFT

31 6025

Residual volume: 1.5 ml

ELS 200 ml (S) syringe

31 6026

Residual volume: 1.5 ml

Packing unit: 50 pieces/carton – 1,500 pieces/pallet



MRS 222 MR-tube system with inline drip chamber (for Scanbag CM-supply bag)

31 7099

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.6 ml
NaCl-Side	36 cm	1.3 ml

Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



Article number

MRS 222 MR-tube system

31 7100

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.6 ml
NaCl-Side	36 cm	1.3 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet



Double head tube system

31 7101

Length from Y-piece: 1.80 cm
 Inner diameter: 2.0 mm
 Total filling volume: 6.9 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet



MR-ELS tube system with valves

31 7104

Inner diameter: 1.7 mm
 Coiled length: 1.80 cm
 Short side: 10 cm
 Long side: 34 cm
 Filling volume: 5.1 ml

Packing unit: 50 pieces/carton



Double head tube system with inline drip chamber (for Scanbag CM-supply bag)

31 7109

Pressure side, ID 2.0 mm	Length	Filling volume
CM-Side	12 cm	0.4 ml
NaCl-Side	36 cm	1.1 ml
Patient line (ID: 1.5 mm)	150 cm	2.6 ml
Suction side, ID 2.7 mm	Length	Filling volume
CM-Side, green	100 cm	5.8 ml
NaCl-Side, white	100 cm	5.8 ml
2 x drip chamber		each one 10 ml

Packing unit: 50 pieces/carton – 2.800 pieces/pallet

CONSUMMABLES : MR

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



MRS 225 MR-tube system

Article number

31 7125

Inner diameter: 2.0 mm

	Length	Filling volume
CM-Side	6 cm	0.2 ml
NaCl-Side	32 cm	1.0 ml

Packing unit: 50 pieces/carton – 4,000 pieces/pallet



MRS 228 MR-tube system, spike

31 7128

1 piece MRS 228 tube system

1 piece SP227 vented spike, luer lock

	ID	Filling volume
Pressure side	2.0 mm	1.5 ml
Suction side	2.7 mm	5.8 ml
Spike	0.4 ml	

Packing unit: 50 pieces/carton – 2,800 pieces/pallet



ES 224/25 patient line with valve

31 8020

Length: 25 cm

Inner diameter: 1.5 mm

Filling volume: 0.4 ml

Packing unit: 100 pieces/carton – 8,000 pieces/pallet



1-way stop cock with 30 cm extension line, luer lock

31 8025

Length: 32 cm

Inner diameter: 2.0 mm

Filling volume: 1.2 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

Double head injectors for magnetic resonance imaging (21 bar/305 psi)

Article number



Patient line 150 cm with terminal check valves

31 8026

Length: 150 cm
Inner diameter: 2.0 mm
Filling volume: 4.8 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Patient line 150 cm with two valves

31 8036

Length: 150 cm + 10 cm
Inner diameter: 2.0 mm
Filling volume: 5.2 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet




ES 224 Patient line

Inner diameter: 1.5 mm

ES 224/80 Length: 80 cm, filling volume: 1.5 ml	31 8081
ES 224/100 Length: 100 cm, filling volume: 1.8 ml	31 8101
ES 224/150 Length: 150 cm, filling volume: 2.7 ml	31 8151
ES 224/200 Length: 200 cm, filling volume: 3.6 ml	31 8199
ES 224/250 Length: 250 cm, filling volume: 4.4 ml	31 8250
ES 224/300 Length: 300 cm, filling volume: 5.3 ml	31 8300

Packing unit: 200 pieces/carton – 11,200 pieces/pallet



DW 229 connector with 3-way stop cock

31 8042

Filling volume: 0.5 ml

Packing unit: 100 pieces/carton

CONSUMMABLES : MR

Double head injectors for magnetic resonance imaging (21 bar/305 psi)

Article number



Patient line with valve

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.6 ml

Length: 200 cm, filling volume: 3.5 ml

Length: 250 cm, filling volume: 4.4 ml

Length: 300 cm, filling volume: 5.3 ml

31 8152

31 8198

31 8251

31 8301

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Patient line with inline valve

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.6 ml

Length: 200 cm, filling volume: 3.5 ml

31 8154

31 8204

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line 150 cm

31 8213

Length: 150 cm

Inner diameter: 1.6 mm

Filling volume: 3.1 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line with 2 valves

Inner diameter: 1.5 mm

Length: 150 cm, filling volume: 2.6 ml

Length: 200 cm, filling volume: 3.5 ml

Length: 250 cm, filling volume: 4.4 ml

Length: 300 cm, filling volume: 5.3 ml

31 8217

31 8218

31 8219

31 8220

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMMABLES : MR

Double head injectors for magnetic resonance imaging (21 bar/305 psi)

Article number



Coiled line with valve

Inner diameter: 1.6 mm

Length: 150 cm, filling volume: 3.1 ml

31 8214

Length: 180 cm, filling volume: 3.7 ml

31 8181

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

31 8252

Length: 300 cm, filling volume: 5.3 ml

31 8302

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line with valve

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

31 8253

Length: 300 cm, filling volume: 5.3 ml

31 8303

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



Coiled line with inline valve

Inner diameter: 1.5 mm

Length: 250 cm, filling volume: 4.4 ml

31 8254

Length: 300 cm, filling volume: 5.3 ml

31 8304

Packing unit: 50 pieces/carton – 5,850 pieces/pallet

CONSUMMABLES : MR

Double head injectors for magnetic resonance imaging (21 bar/305 psi)



Coiled line 235 cm with Y-connector
with check valve

Article number

31 8235

Length: 235 cm + 1 x 10 cm + 1 x 12 cm
Inner diameter: 1.6 mm
Total Filling volume: 5.6 ml

Packing unit: 50 pieces/carton – 5,850 pieces/pallet



DW 226 3-way stop cock

31 9085

Filling volume: 0.3 ml

Packing unit: 300 pieces/carton



SP 227 vented spike, luer lock

31 9095

Filling volume: 0.4 ml

Packing unit: 50 pieces/carton



Minispike, green, luer lock

31 9099

Filling volume: 0.2 ml

Packing unit: 300 pieces/carton